



Horizon® - Specifically designed for mobile integration

Simplify Imaging

- **Print Versatility-** Horizon's unique ability to print film, color and grayscale paper in the same compact imager provides greater versatility and savings over "film-only" imagers. Horizon provides an overall total lower cost solution on film, maintenance and service than any competing imager, making it the ideal solution for any mobile application.

Rugged Design

- **Reliability-** Codonics adapted rigorous testing methods used in the aerospace industry, ran Horizon through millions of print cycles, as well as crash, impact and drop testing. Horizon withstands the significant, prolonged exposure to vibration and rapid transitions in temperature with no degradation in image quality.

Improve Workflow

- **Easy to Site and Service -** Horizon's ultra-compact footprint allows it to be placed on a desk, counter or in any mobile application. It weighs less than 70 lbs. and features standard power requirements with plug-and-play simplicity.
- **Speed -** With the fastest time-to-first-print, Horizon prints up to 100 films an hour.





Horizon®

Specifically Designed for Mobile Applications.

Codonics engineered Horizon for the most severe mobile conditions with test data to back it up. Horizon is rugged enough to withstand being transported into some of the harshest places on earth. Its durable cabinet mounts in any mobile coach while the docking tray secures the imager and provides a quick release should a service event occur. Horizon's unique swap service and smart card technology eliminates the need for on-site service calls. No other imager offers this seamless integration in a mobile environment.

Horizon is easily configured to accommodate multiple sites on one route using our exclusive built-in Zip® Drive to download configuration data specific to each site. There is no limitation to the number of configurations that can be created.

Optional docking tray and cabinet with storage




Mobile Applications

- Direct or Network Connection
- High Throughput
- Compact and Easy to Site
- Reduce Costs
- Video and PostScript Compatible
- Satisfy Demands from Referring Physicians

Referring Physicians

- Clinical Review
- Patient Consults
- Patient Medical Files
- Surgical Planning
- Patient Marketing



17991 Englewood Drive
 Middleburg Heights, OH 44130 USA
 (440) 243-1198
 (440) 243-1334 Fax
 Email info@codonics.com
 www.codonics.com

For more information contact Codonics at 800-444-1198 or 440-243-1198 and ask for our Mobile Sales Manager. Copyright © 2005 Codonics, Inc.