

Safe Label System®

Administration Tool (AT)

EN - English

Overview

Never before has it been easier to empower every medication preparation area with pharmacy-level oversight and compliance. Codonics Safe Label System (SLS) Point of Care Station (PCS) sits wherever medications are prepared and utilizes the integrated formulary database to increase patient safety and labeling compliance. This database, created using the SLS Administration Tool (AT), is the system's management software that allows pharmacy to easily create and control a drug formulary, including site-approved medications, concentrations and rules for use, expanding pharmacy's scope of influence to the point of medication preparation. The AT enables a hospital pharmacy to:

Features/Benefits

- Quickly and easily create a formulary that incorporates hospital-approved drugs, merged with hospital policies, together with the Joint Commission (TJC) requirements, and American Society of Anesthesiology (ASA) and ISO labeling guidelines
- Control deployment of the formulary, configuration, SLS user's packages and software updates to any SLS PCS
- Load site-approved, customized formularies on each SLS PCS so that they are available to clinicians at the point of preparation and administration to safely verify, identify and confirm medications to increase patient safety and reduce medication errors
- Reduce the amount of time spent on manual paperwork maintaining 340B compliance when used in conjunction with SLS-WAVE SLS includes a one-year subscription to Lexicomp™ Basic Drug Database. To populate the MDD with drug records users can import the Lexicomp database, their own database or add drugs manually.

Specifications

Recommended - Requirements for AT Server

Hardware: Intel Quad Core (i5, i7 or Xeon) with 4 cores, 8 GB RAM

Hard Drive Space: Minimum 5 GB free space for application. Recommended 30 GB total at install

Computer OS: Windows® 7, Windows 10, Windows Server 2008, 2012R2 and 2016. 64-bit required.

Recommended - Requirements for AT Client

Browser Support: Internet Explorer IE 8 or newer

Screen Resolution: Minimum 1440 x 900; use default browser fonts

Computer OS: Windows® 7 or newer

Applications: Excel® 2003 or greater to open reports

Single Sign On: AT User Management using Microsoft® Active Directory