

Codonics Horizon Operating Software

Technical Brief

Version 3.0.2 Release Notes

Summary

Horizon® 3.0.2 software now includes the US FDA Unique Device Identifier Global Trade Item Number (UDI GTIN). UDI will allow industry and regulatory authorities to more rapidly identify medical devices.

Products Affected

Imagers with a serial number prefix of 53C, 54C, or 55C are affected. Please contact your Codonics representative for further information.

The Imager Software Version can be determined using the Control Panel by viewing the Imager Info screen from the Utilities menu.

NOTE: Horizon Imagers with serial number prefix 53C, 54C, or 55C cannot be downgraded to any previous software version below 3.0.0. Imagers with a serial number prefix 53C, 54C, or 55C cannot have their software upgraded to 4.0.0 or higher.

New Features (Adaptive Updates)

- ◆ The FDA UDI GTIN (14-digit number) is now included in the Horizon software. It can be found in the Quick Info Screen.

Known Common Issues

- ◆ Previous issues can be found in previous Release Notes.

Known Uncommon Issues

- ◆ Previous issues can be found in previous Release Notes.

Technical Support

If problems occur during software installation, please contact Codonics Technical Support at any time.

Phone: +1.440.243.1198

Email: support@codonics.com

Website: www.codonics.com

Get it all with just one call
800.444.1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.
Copyright © 2016 Codonics, Inc. Printed in the U.S.A.
Part No. 901-254-003.01.



17991 Englewood Drive
Middleburg Heights, OH
44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com

Codonics Trading Co, Ltd.
317 Xianxia Rd. Building B
Unit 1412
Changning Dist., Shanghai
P.R. China, 200051
86-21-62787701
86-21-62787719 Fax

Codonics Limited KK
AQUACITY 9F,
4-16-23, Shibaura
Minato-ku, Tokyo,
108-0023 JAPAN
81-3-5730-2297
81-3-5730-2295 Fax