



CODONICS

We bring the future into focus

Infinity™

Medical Image Server

IHE Integration Statement

Part Number: 900-457-001

Document / Software Version: 1.3.0

Date: 5/11/2009 3:38:00 PM

Revision	Date	Author	Description of Changes
01	1 April 2009	Rich Edwards	Initial Release
02	5/11/09	Ross Goodman	Updated copyright date and version date.

IHE Integration Statement			
Vendor	Product Name	Version	Date
Codonics Inc.	Infinity Medical Image Server	V1.3.0 and later	3 April 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options	Implemented
Consistent Presentation of Images	Image Manager/Image Archive	N/A	
Access to Radiology Information	Image Manager/Image Archive	N/A	
Key Image Note	Image Manager/Image Archive	N/A	
Simple Image and Numeric Report	Report Repository	N/A	
Reporting Workflow	Image Manager/Image Archive	N/A	
Evidence Documents	Image Manager/Image Archive	N/A	
NM Image	Image Manager/Image Archive	N/A	
Mammography Image	Image Manager/Image Archive	N/A	
Audit Trail and Node Authentication	Secure Application	Radiology Audit Trail	
	Audit Record Repository		
Internet address for vendor's IHE information: www.codonics.com/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.codonics.com/Support/Infinity/pdf/Infinity_DICOM_Conformance_Statement.pdf		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	